

The Iron-Rhodium System

Modification of the Magnetic Phase



Member of the Helmholtz Association

Nico Spahn | Institute for Ion Beam Physics and Materials Research – Interface Magnetism | www.hzdr.de

Outline


1. Motivation - The FeRh binary alloy
2. Introduction
3. The new measurement system
4. Outlook



Copyright by Time Warner Entertainment

Motivation: Heat assisted magnetic recording (HAMR)

- *Higher storage Density
- *Lower writing speed




Source Seagate: Whoever said disk drive technology and supply was settling into boredom? No way, Jose!

Outline

1. Motivation
2. **Introduction**
3. The new measurement system
4. Outlook

Introduction: The Iron-Rhodium System




Phasediagram of the FeRh-Alloy.

Swartzenhuber, L.J., Bulletin of Alloy Phase Diagrams, 1984


Near equiatomic composition:

- B2-structure (bcc based, AFM, FM)
- A1-structure (fcc based, PM)
- Thermally induced phase transition




Introduction - The Iron-Rhodium system

AFM ground state



- α'' -phase (B2)
- G-type AFM order
(magnetic monolayers along 111 plane \rightarrow GMR)
- $\text{Fe} \approx 3,3 \mu_B$
- $\text{Rh} = 0 \mu_B$

FM high temperature state




- α' -phase (B2)
- FM order
- $\text{Fe} \approx 3,2 \mu_B$
- $\text{Rh} = 0,9 \mu_B$

Introduction - Basics for the measurement


Resistance change at Curie point

- Occurs in ferromagnets
- Electrical resistance below Curie-temperature smaller than in non-ferromagnets
- Above the Curie-Temperature ferromagnets become paramagnetic and behave like normal conductors



Comparison of the electrical resistance of the non-ferromagnetic Pd and the ferromagnetic Ni.
(R.Gross, WMI Lecture Spintronics , 2005)

Introduction - Basics for the measurement



Magnetic phase transition and negative magnetoresistive effect of a 50 nm thick FeRh-Layer.

Suzuki, I., et al., Journal of Applied Physics, 2011

AFM-FM phase transition

- around 350 K

Phase transition detection

- Resistance
- Magnetization
- Giant Magnetoresistance



Between adjoined Fe layers

$\text{Fe}\uparrow\text{-Rh-Fe}\downarrow - R \uparrow$

$\text{Fe}\uparrow\text{-Rh-Fe} \uparrow - R \downarrow$


Introduction - Ion beam induced magnetization changes

Ion irradiation of 80 nm B2 FeRh films



Ion irradiation

- Change in magnetization as a function of fluence
- Critical value of fluence exists
- Above critical value change of lattice structure from B2 to A1



N. Fujita, Journal of Applied Physics, 2010

Outline

1. Motivation
2. Introduction
3. The new measurement system
4. Outlook

The New Measurement System

How ???

FeRh properties

- AFM-FM transition due to temperature at 350 K
- AFM-FM transition due to ion irradiation
- Detection by resistance change with temperature
- Detection by Giant Magnetoresistance measured as function of external field at $T > 350\text{K}$




In-situ system with

- Temperature variation
- Resistance measurement
- GMR measurement (external field required)
- Ion irradiation and measurement of the fluence



The New Measurement System

Material Selection





The New Measurement System

Construction



- Heating up to 400+ K
- Chamber reaches high vacuum conditions
(leakage rate $10^{-10} \text{ mbar} \cdot \text{l/s}$)
⇒ important for running the ion gun

Construction of the Measurement System



Construction of the Measurment System

Safety - High voltage is not your friend !!!



Ion irradiation

- Creates charging effects of the sample stage about 5 kV
- No transfer of the charge to any touchable parts
- Secured grounding of this charge



PEEK as high performance polymer with a dielectric strenght of **75 kv/mm**

Conclusion

Possible measurements

- Measuring the GMR via layer resistivity with applied external magnetic field
- Measuring the phase transition temperature vs. ion flux via the layer resistivity
- Under high vacuum conditions

Acknowledgements



You don't reach your goals without
some friends.

Dr. Kay Potzger
Andreas Henschke
Alireza Heidarian

Thank you !