


WP 3

“Time-resolved IR spectroscopy on proteins”




PIDID Meeting
December 16th 2011

Dr. Christian Bauer


Guideline

1. Motivation for THz spectroscopy of proteins
2. WP3 "Time-resolved IR spectroscopy on proteins"
 - Advanced pump/probe experiment
 - Model protein Bacteriorhodopsin: structure, function and spectroscopic properties
 - Experimental setup of advanced pump/probe experiment
3. Results: FIR difference spectroscopy on Bacteriorhodopsin (preliminary work)

Spectroscopy on Proteins in the Infrared



Spectroscopy on Proteins in the Infrared




mid infrared: standard method

Detection of change of vibration of

- Chromophores
- Protein backbone
- single amino acid side chains

→ structural analysis on molecular level

Spectroscopy on Proteins in the Infrared



far infrared (THz):
new method of light induced difference
absorption spectroscopy

Detection of change of vibrations of

- Chromophores
 - Protein backbone
 - H-bonding network
 - Low energy torsional modes
- structure analysis on molecular level

mid infrared: standard method


Detection of change of vibration of

- Chromophores
 - Protein backbone
 - single amino acid side chains
- structural analysis on molecular level

Motivation for Spectroscopy on Proteins in the Far Infrared

- Assignment of far IR modes helps to understand the functionality of proteins
- In the far infrared collective vibrations of the chromophore and the H-bonding network are excited
- Far Infrared radiation can be used to perform special multi-wavelengths pump/probe experiments

Classic Pump/Probe Experiment



- **pump pulse**

stimulates a process,
reaction, phase transition...


-> preferably by photons!

- **probe pulse**

measures change of
materials properties at
different times

Modified Pump/Probe Experiment: Control of Photoreaction

Pump with 2nd
wavelength
(FEL)



- **pump pulse**
stimulates a process,
reaction, phase transition...


-> preferably by photons!

- **pump pulse with 2nd
wavelength**
changes the photoreaction

-> control of reaction


- **probe pulse**
measures change of
materials properties at
different times

Setup of modified Pump/Probe Experiment



two pump beams: VIS and THz-FEL
 one probe beam: IR (FTIR-Globar)
 time-resolved stepscan technique

Time-resolved FTIR measurement: "Stepscan" technique





- Preliminary stepscan experiment without sample
- Proof of principle with "chopper experiment" according to manual

Attenuated Total Reflection (ATR) technique


Window material in transmission setup not transparent in the FIR
(KBr below 25 μm , CaF_2 below 12 μm , BaF_2 below 15 μm transparent)

- + Setup vacuum-tight
- + Sample preparation possible as film
- Lower signal-to-noise ratio as in transmission




Modified Pump/Probe Experiment Biomolecule: Bacteriorhodopsin (bR)

3D structure of Bacteriorhodopsin




Photocycle of Bacteriorhodopsin

simplified photocycle of Bacteriorhodopsin





M-BR difference absorption spectrum




Photocycle of Bacteriorhodopsin

simplified photocycle of Bacteriorhodopsin




Difference Absorption Spectrum of a Bacteriorhodopsin Film in the Mid Infrared



- Characterization of accumulated intermediate using spectra from literature: M state
- Optimization of sample preparation


Light induced Difference Absorption Spectrum of a Bacteriorhodopsin Film in the Far Infrared



Difference bands in the far IR


Ratio of the amplitude of difference bands in the mid IR and the far IR $\approx 10:1$
 \rightarrow experiment in the FIR takes ≈ 100 times longer than in the mid IR

Light-induced Difference Absorption Spectrum of a Bacteriorhodopsin film in the Far Infrared




- positive and negative difference bands
- without sample: no difference bands, stronger noise at 380 cm^{-1}


Tentative Assignment of the Negative Difference Bands



all-trans Retinal in bR,
protonated Schiff base



all-trans Retinal without apoprotein



→ Similarity in band pattern

Samples comparable despite differences in:




- Structure
- Chromophore environment

Tentative Assignment of the Negative Difference Bands

Difference bands of bR film		all- <i>trans</i> Retinal in PE	
positive	negative	this work	Gervasio <i>et al.</i>
482 cm ⁻¹			
	468 cm ⁻¹	487 cm ⁻¹	487 cm ⁻¹
451 cm ⁻¹			
	441 cm ⁻¹	467 cm ⁻¹	466 cm ⁻¹
433 cm ⁻¹			
	424 cm ⁻¹	409 cm ⁻¹	409 cm ⁻¹
408 cm ⁻¹			
	386 cm ⁻¹	377/370 cm ⁻¹	377/369 cm ⁻¹
366 cm ⁻¹			
	327 cm ⁻¹	343 cm ⁻¹	331 cm ⁻¹
315 cm ⁻¹			
	283 cm ⁻¹	300/290 cm ⁻¹	301 cm ⁻¹
258 cm ⁻¹			
	216 cm ⁻¹	179 cm ⁻¹	180 cm ⁻¹
170 cm ⁻¹			
	161 cm ⁻¹	163 cm ⁻¹	162 cm ⁻¹

Tentative Assignment of the Negative Difference Bands

Negative difference bands of bR film	all- <i>trans</i> Retinal in PE	Assignment of bands by calculations Gervasio <i>et al.</i>
468 cm ⁻¹	487 cm ⁻¹	C8C9C10 bending + C8C9C10-C19 bending
441 cm ⁻¹	467 cm ⁻¹	
424 cm ⁻¹	409 cm ⁻¹	C8C9C19 bending + C12C13C20 bending
386 cm ⁻¹	377/370 cm ⁻¹	
327 cm ⁻¹	343 cm ⁻¹	
283 cm ⁻¹	300/290 cm ⁻¹	
216 cm ⁻¹	179 cm ⁻¹	
161 cm ⁻¹	163 cm ⁻¹	C13-C20H₃ torsional



Conclusion

- Development of an ATR-Setup for spectroscopy in the far IR
- Measurement of light-induced difference spectra of a Bacteriorhodopsin film in the far IR from 490 to 150 cm^{-1}
- Negative difference bands tentatively assigned to vibrational modes of all-*trans* retinal

Outlook

Band assignment useful knowledge to perform two-colour-pump/probe experiments:

Re-isomerisation of Bacteriorhodopsin's retinal by resonant excitation in the far IR possible → **control of photocycle**

Prokhorenko, V. I. *et al.* Science **313**, 1257 (2006)

„coherent control“ using fs pulses in the NIR: control of isomerisation yield of Bacteriorhodopsin shown

Acknowledgement



University Bielefeld / FU Berlin

Prof. Dr. J. Heberle

Dr. T. Kottke

Dr. S. Kakorin

Dr. K. Ataka

C. Pudel

I. Ehring



HZDR

Dr. M. Gensch

Dr. W. Seidel

Project funded by the BMBF „PIDID“ # 05K10KEB



Thank you for your attention